Eric Pitman Summer Workshop in Computational Science




Jeanette Sperhac & Amanda Ruby



Introducing the Workshop Project


Here's what we'll cover:


- The story of the HWI protein crystallization data
- The Questions
- So what's a classifier?
- Inside the dataset
- The Project in RStudio
- Exploring the Proteins
- What you'll need



Proteins

Proteins are large biological molecules, composed of long chains of amino acids






These chains fold into complex structures, allowing the protein to perform biological tasks

STRUCTURE AND FUNCTION ARE VERY CLOSELY RELATED


X-Ray Crystallography

- Protein crystals are bombarded with x-rays. The x-rays are diffracted, giving structural biologists a way to determine structure
- Herbert A. Hauptman developed the mathematical model to convert reflections to molecular structures
- The process requires *protein crystals*



Crystallization


- Trying to get a protein to "crash out" of solution, in an orderly fashion
- Requires a precise set of conditions, which vary from protein to protein
- Included in this precise set of conditions: a chemical cocktail
- Cocktails are combined into screens (==generations)




Introducing the Workshop Project

To achieve this (protein structure):


We must get this (pure, crystallized protein):






Protein Data Bank ID: 3dm3

HWI Protein Crystallization Data






Which sample contains a crystallized protein? (Can we make this decision for 3 million samples?)

Real Data, Local Data



HWI lab work





CCR processing

1. Laboratory Work: HWI



- A protein is placed in 1536 wells on an experimental plate.
- A different chemical cocktail is added to each well.
- Photos of each well are taken at different timepoints.

Which proteins crystallized in which cocktails, at which timepoints?


2. Automatic Classifier: CCR



An automatic classifier helps identify crystallized samples:

- Compute measures that describe the photo of each sample.
- Classifier assigns probability that each sample contains crystallized protein.

3. But...




The classifier is only ~70% accurate, so:


Human expert classifies each sample. This is the final word.

39936 Experimental Records

We have a photograph of each experiment:

- 13 proteins
- In 1536 chemical cocktails
- At 2 different time points







39936 Experimental Records

We have information about each one:

- Classifier score (value between 0 and 1)
- Human assignment (truth): crystal/non-crystal
- Some other stuff...





Automated Classifiers

From: cheapsales@buystufffromme.com

To: ang@cs.stanford.edu

Subject: Buy now!

Deal of the week! Buy now!

Rolex w4tchs - \$100

Medicine (any kind) - \$50 Also low cost M0rgages

Also low cost Murgages

available.

Class 1 Span

From: Alfred Ng

To: ang@cs.stanford.edu Subject: Christmas dates?

Hey Andrew,

Was talking to Mom about plans for Xmas. When do you get off

work. Meet Dec 22?


Alf

Non-spom

Class 2

This classifier has assigned data points into two classes, 1 and 2.


Automated Classifiers



How does this work?

- Classifier is run multiple times on each data point
- The results are converted to a probability that a crystal was seen

Classifier Performance




Our classifier assigns scores that lie between 0 and 1.

We can overplot two kernel density plots:

- One of human-classified crystal examples
- One of human-classified non-crystal examples

...both with classifier score on the x axis

Classifier Performance




Pick a threshold (cut point) for classifier scores:

- above it, consider all examples to be positive (crystal)
- below it, consider all examples to be negative (not crystal).

No matter where we cut, we get some classifications wrong!


Classifier Performance



After we pick a threshold, we can see four types of outcomes:


- True and False Positive (crystal)
- True and False Negative (not crystal).

ROC Plots and Classifiers



- Evaluate the performance of a binary classifier
- Vary the threshold or cut point, vary the sensitivity/specificity ratio

Sensitivity vs. Specificity



- Specificity (true negative rate): the ability to classify noncrystals as noncrystals.
- Sensitivity (true positive rate): the ability to classify crystals as crystals

Meet Your HWI Dataset

The data frames:

experiment: 39936 rows

• sample: 13 proteins

drop: 3930 rows (1536 chemical cocktails)

expUrl: 39936 rows (one per experiment)

 $13 \times 1536 \times 2 = 39936$ proteins cocktails timepoints experiments

Experiment Data: 39936 Rows

experiment data frame describes the crystallization experiments:

- read_no: 39936 values, identifies experiment
- sample_id, plate_no: 13 values, identifies protein
- week_no: 2 timepoints
- cocktail_no: 1536 unique values
- human_crystal: crystal or not?
- class3_crystal: probability of crystal

```
13 \times 1536 \times 2 = 39936 proteins cocktails timepoints experiments
```

Sample Data: 13 Different Proteins

sample data frame describes the proteins:

- Identify each protein:
 - sample_id, p_number, targetdb_status, name
- experimental_molecular_weight
- concentration of the protein
- seq: sequence of amino acids in protein
- seq_len: number of amino acids
- targetdb_ref: PDB (Protein Data Bank, www.rcsb.org)
- targetdb_status: unique protein code on PDB

Drop Data: 3930 Rows


drop data frame describes the components of the chemical cocktails:

- Identify each chemical cocktail (1536 in total):
 - cocktail_no + solution_component_no
- concentration of the solution component
- name of the solution component
- ph: pH of the solution component

expUrl Data: 39936 Rows

expUrl data frame contents:

- read_no: identifies experiment
- image_url: URL of the experiment image





HWI data are text files with csv format.

- Formatted to load straight into an R data frame
- Columns are labeled by name

csv == comma separated values




RStudio and csv: Two ways to load data

- 1. RStudio Workspace:
- Select Import Dataset: From Text File
- Select a .csv file to Open
- Use Heading=Yes

- 2. Or, from the command line:
 - Set the Working Directory
 - Load command:

```
> drop=read.csv("drop.csv")
```

Exploring Proteins



Protein Data Bank ID: 3dm3

- BLAST algorithm
 - Protein sequence comparison
 - R package: bio3d
 - NCBI BLAST web interface
- Protein Data Bank

Exploring Proteins: bio3d BLAST


blast.pdb() returns:

- bitscore
- evalue
- mlog.evalue
- pdb.id


bio3d package documentation: http://thegrantlab.org/bio3d/html/index.html

plot.blast() returns:


Summary and plot of BLAST hit statistics




Exploring Proteins: NCBI BLAST



Exploring Proteins: rcsb.org



Exploring Proteins




Check the data frame expUrl for links to the experimental images.


Can you identify the crystallization state?

| | read_no | image_url |
|---|----------------------------|---|
| 1 | X0000095890696200801181226 | http://xtuition.ccr.buffalo.edu/image-data/X000009589/X000009589200801181215/X0000095890696200801181226.jpg |
| 2 | X0000095890384200801181231 | tuition.ccr.buffalo.edu/image-data/X000009589/X000009589200801181215/X0000095890384200801181231.jpg |
| 3 | X0000095890001200801181235 | ccr.buffalo.edu/image-data/X000009589/X000009589200801181215/X0000095890001200801181235.jpg |
| 4 | X0000095890488200801181228 | ion.ccr.buffalo.edu/image-data/X000009589/X000009589200801181215/X0000095890488200801181228.jpg |
| 5 | X0000095890504200801181228 | ion.ccr.buffalo.edu/image-data/X000009589/X000009589200801181215/X0000095890504200801181228.jpg |
| 6 | X0000095890926200801181222 | tion.ccr.buffalo.edu/image-data/X000009589/X000009589200801181215/X0000095890926200801181222.jpg |
| | | |

The questions

- 1. How does the automated *classifier* perform?
- 2. What trends do we see in *timeseries* data?
- 3. What trends do we see in crystallization of the proteins across the different *cocktails*?
- 4. What can we learn about the *proteins*?





Project: you'll need...

- Data load
 - read.csv() R function
- File transfer: WebDAV
 - Copy files to your workstation from hpc2.org
 - Useful for creating your presentation slides
- factorify() function
 - Apply to a data frame after subsetting; updates factor levels
 - source generalFactorifyFunction.R
- Optional: Jmol protein viewer on hpc2.org

1. Classifier Performance Focus

- Data file: experiment.csv
- Epi R package: ROC()
- R functions:
 - sapply()
 - Kernel density estimation: density()
 - data.frame()
- R graphics:
 - hist(), plot()
- Provided function compareDensityPlots()

2. Timeseries Crystallization Focus

- Data files:
 - experiment.csv
 - sample.csv
- R functions:
 - merge(), length(), table(), round()
- R graphics:
 - boxplot(), hist(), pie(), barplot(), legend()

3. Cocktail by Protein Focus

- Data files:
 - Experiment.csv
 - Sample.csv
 - Drop.csv
- R functions:
 - table(), unique(), sum(), length(), which(), dim(), merge()
- R graphics: boxplot(), pie(), barplot()

4. Protein Exploration Focus

- Data files:
 - expUrl.csv, experiment.csv, sample.csv
- PDB (protein data bank): www.rcsb.org/pdb
- NCBI BLAST web interface
- R "bio3d" package: blast.pdb(), plot.blast()
- R functions:
 - plot() for scatterplots
 - Im() linear model to fit data to a line
 - abline() plot a linear model
 - cor() examine correlations